



Smitech World's 1st Covid-19 Protective System to Mitigate Risks of Infection



- Viruses are protected by a layer of lipid (fat) and protein membrane.
- Thrive under low temperatures.
- Infects by attaching its protein to ACE2 (human cell) and release genetic material of virus.
- WHO confirmed coronavirus is spread by droplets and airborne pathogens from infected people.
- Asymptomatic cases have low viral dose.


Research Background

- Virus is a self-assembled nanoparticle in which the weakest link is the lipid (fatty) bilayer.
- Studies have shown that the lipid bilayer of virus will melt at 56°C after 15 minutes of exposure (Source: WHO)
- Pure oxygen will boost the cell immune system and assist in breathing difficulty for mild cases.
- Heated O₂ breaks the lipid (fat) surrounding viruses, causing them to “melt” and perish.



Future Outlook of Coronavirus Pandemic

1. WHO has confirmed that many people will not have the old normal for long time to come
2. Predicted 3 Waves similar to Spanish Flu lasting 2 years. The second wave has the highest fatalities
3. CNN / CDC confirmed 40% of Covid-19 tested positive are from asymptomatic cases
4. Inhaling heated humidified oxygen can reduce viral load by 80%. Many successful reports using heat



Spanish Flu death toll across three waves
2nd wave being the highest fatalities.

Protect your Family & yourself with our System



Health is Wealth

SMITECH (ASIA) PTE LTD
 75 BUKIT TIMAH ROAD, #06-13 BOON SIEW BUILDING
 SINGAPORE 229833 TEL: +65 6339 5277 FAX: +65 6339 0253
 E: sales@smitechasia.com.sg or smitech@pacific.net.sg
 W: www.smitechasia.com.sg



Your Partner
 For Healthy Lifestyle

Covid-19 Protective System Using Warm Oxygen Inhalation Treatment



How our system works

1. Reduce viral load by 80% after each treatment cycle
2. Protect dentists and patients against infection
3. Prevent Covid-19 at home, clinics, isolation centres and hospitals
4. Assist in breathing difficulty for high risk patients
5. Safe to inhale warm humidified oxygen below sauna & steam bath temperature and duration

Treatment Procedure

- Inhale pure heated humidified oxygen (60°C-69°C) for 15 to 20 minutes then rest for 40 minutes. This is one cycle.
- This operating temperature is far below that of saunas and steam baths. Recommended 6-8 cycles each day to neutralise all virus load.
- Pure oxygen will provide patient's body with the fuel it needs to build new cells and fuel one's immune system.
- Monitor heart rate and SpO2 level using a pulse oximeter during treatment process.

This system consists of :

- Smitech oxygen generator
 - External oven with heating coil
 - Temperature controller with mask sensor
 - 10 pcs reusable oxygen mask
- Optional:
- Storage rack
 - Fingertip pulse oximeter
 - Electronic timer
 - Hydrogen water generator to remove harmful free radicals and provide anti-inflammatory benefits

Benefits of our System

- Oxygen is needed to construct 7 hundred billion replacement cells each day and it is essential for our immune system.
- Our immune system kills pathogens and "fuels the cells that make up our body's defences against viruses and other invaders."
- This non-invasive treatment can shorten 14-day quarantine period by 10 days if tested negative.
- It can disinfect patients before dental treatment & protect dentists & staff
- It can disinfect all high risks people in hospitals, homes, sports, crowded places
- It has a very silent operation.
- Small foot print of 50cm x 80 cm space.

Options

- **Rental : On request (minimum of 2 months)**
- **Sales : On request**
- **Warranty: 2 years or 10,000 working hours**

Operation: AC230V 50Hz 1200W
CE approved / ISO 9001:2015



Health is Wealth

SMITECH (ASIA) PTE LTD
75 BUKIT TIMAH ROAD, #06-13 BOON SIEW BUILDING
SINGAPORE 229833 TEL: +65 6339 5277 FAX: +65 6339 0253
E: sales@smitechasia.com.sg or smitech@pacific.net.sg
W: www.smitechasia.com.sg



Your Partner
For Healthy Lifestyle