

# Portable Oxygen Concentrator with battery



**OLV-C1 portable oxygen concentrator with *Battery***  
is the most lightest, smallest and quietest oxygen  
concentrator for travel and ambulatory oxygen needs.

#### DISTRIBUTOR:

**SMITECH (ASIA) PTE LTD**

75 BUKIT TIMAH ROAD, #06-13 BOON SIEW BUILDING

SINGAPORE 229833 TEL: +65 6339 5277 FAX: +65 6339 0253

E: [sales@smitechasia.com.sg](mailto:sales@smitechasia.com.sg) or [smitech@pacific.net.sg](mailto:smitech@pacific.net.sg)

W: [www.smitechasia.com.sg](http://www.smitechasia.com.sg)



# FEATURES

1. Continuous flow, 1-5L/min adjustable
2. Light weight: 5kg
3. Three power supplies: 220V/110V, 12V(Car use), Rechargeable Lithium-ion Battery
4. Come with carrying bag and portable trolley. It is more convenient to use at outdoor
5. Low noise, less than 40db(A)
6. Anion function



## SPECIFICATIONS

<b>Size</b>	310*165*310mm
<b>Weight</b>	machine 5kg, battery 1kg
<b>Flow Settings</b>	1-5LPM
<b>Purity Oxygen</b>	≥90%±3%(1LPM)
<b>Output Pressure</b>	40Kpa-60Kpa
<b>Sound Level</b>	less than 40db(A)
<b>Power</b>	1. AC Power 220V±22V, 50Hz±1Hz/110V±15V, 60Hz±1Hz 2. DC Power 12V 3. Rechargeable Battery
<b>Power Consumption</b>	90W
<b>Battery Duration</b>	60mins
<b>Battery Recharge Time</b>	4 hours recharge time to achieve 90% capacity